

**A flexible Cardiovascular Laser-Catheter
INSTRUCTION FOR USE (IFU)
GENERAL DESCRIPTION**

The **HypertenoLas®** is a tripolar (8F) laser catheter designed for percutaneous transluminal renal and pulmonary denervation (RD and PD). It can be used for:

- intravascular catheter directed **laser application**
- **selective** anatomically or **mapping guided** retrocardiac, and perivascular nerve modulation / denervation
- Low and high frequency stimulation (LFS, HFS)
- infusion of liquids

For its use the following additional devices are needed:

- an 8.5 French steerable sheath e.g., AGILIS,
- the laser **CardioVascLas®**, 1064, LPS,
- the Rolling pump **IriFlowLas®**
- a low and high frequency current stimulator

Laser catheter flexibility allows for a safe and quick catheterization of the right and left renal artery from the abdominal aorta or of the pulmonary artery from the right ventricular outflow tract.

Laser application is performed under **normothermic conditions**. During laser application the **catheter itself is not heated up**. Laser induced absorptive heating is tissue specific and selective. Sympathetic nerves are **selectively heated up** and are coagulated, with minimal and reversible thermal damage to the vessel walls when using an appropriate energy setting. Light absorption of the transparent vessel intima and media are minimal.

Protection of the optical fiber in the catheter lumen and continuous saline irrigation allows for a **non-contact** mode of laser application. Saline irrigation creates a clear medium for the laser light.

Saline irrigation at room temperature of 18°C has a cooling effect on the inner vessel wall. Whereas deep penetration of the laser light will achieve an enough high level of temperature, needed for thermal damage to the perivascular sympathetic nerves located in a circular network around the vessels.



Perivascular laser denervation is **not thrombogenic**. Laser application under normothermic conditions and continuous saline irrigation avoids intravascular and mural thrombus formation. There is no intravascular temperature increase.

The **divergent laser beam** produces circular radiation of a vessel segment of 4-5 mm. The sympathetic nerves around the vessels are heated up selectively and are coagulated. The length and depth of laser lesion hinders regeneration of permanently inactivated nerves.


The **three 5 mm long cable electrodes** mounted at **interelectrode distances 2mm** longitudinally on the catheter head allow for **LFS** and for **HFS-mapping** guided stimulation and so for localization of perivascular sympathetic and vagal pulmonary innervation.

Online monitoring of pulmonary electrical potentials is performed **during LFS and HFS** without interfering with the electrophysiologic monitoring principles, **without hum** in the ECGs also **during laser application**.


Laser application can be performed without electrical hum in the intracardiac electrograms during LFS / HFS.



Renal nerves originate from proximal ganglia, converge on the renal artery, are closer to the artery distally: **red circles**, where laser application should be aimed at. Accessory arteries are common and are always innervated. RD requires complete treatment of **distal main renal artery, and accessories**.



ECG recordings during intra-pulmonary artery (PA) Low Frequency (LFS) and High Frequency Stimulation (HFS) mapping showing:
A LFS with a cycle length of 333 ms (180 bpm) at 10V and 1 ms in the right PA results in atrial capture indicated by the white vertical arrows.
B No atrial capture during LFS at another point in the right PA
C HFS at the same point as in B with 33Hz (10V, 1ms) for 10s results in a shortening of atrial cycle length from 430 ms to 400 ms.
D After laser application aimed at this point there is **no change** in the atrial cycle length during HFS. RAE = Right Atrial Electrogram



Right: The laser catheter is placed through the steerable sheath (lower arrow) with its tip (upper arrow) orientated towards the pulmonary artery trunk

APPLICATIONS

Indications

The **HypertenoLas®** is designed for renal and pulmonary artery denervation by applying 1064nm laser light, for attenuation/modulation of the nerve overactivity, mainly to reduce high blood pressure. Preferably to **normalize** the systemic or pulmonary blood pressure permanently.

Via its conically shaped optical fiber the vessel walls can be radiated in a circular fashion with laser light. The divergent laser beam allows for circular radiation at lower power density and so helps avoid thermal damages of the vessel wall.

For the RSD or PAD procedures an aortogram or pulmonary artery angiogram must be performed, and for RSD the anatomical dimensions of targeted vessel segments must be calculated. Radiation should be **limited to 15W**, and application times adapted to the irradiated surface of the vessel inner lumen to achieve a safe and effective perivascular denervation. Examples:

For RSD 15W/15s (2-4x15 s) and for PAD 15W/15s (225J), and repeated until unchanged atrial cycle lengths during atrial catheter stimulation are monitored (see page 4).

Contraindications

The use of the **HypertenoLas®** may be contraindicated if there is a known or suspected obstruction in the vessel access or a vessel spasm. Acute and severe chronic diseases, especially a substantially reduced left ventricular function, obstruction of coronary arteries and angina pectoris, increase the procedure risk and may contraindicate the use of the **HypertenoLas®**.

SAFETY NOTES

Reuse

The **HypertenoLas®** is designed for single use only. Cleaning or reuse can result in serious complications. **LasCor** will not be responsible for any direct or consequential damages or expenses which result from cleaning or reuse.

Sterilization

The catheter **is** shipped after sterilization with ethylene oxide (EO) gas. Do not use products from opened or damaged packaging. Under appropriate storage conditions, we guarantee sterility in undamaged packaging until the expiry date (use before date).

External Interference

There is no external interference with the laser catheter **HypertenoLas®**. Sterile products should be stored at humidity of 45-70% and at temperatures of 18-25°C. They should not be exposed to direct sunlight and must be used before the expiry date on the packaging.

Side Effects and Complications

Despite correct handling of the **HypertenoLas®** complications may occur such as:

- Vessel spasm and local bleeding with subcutaneous Hematoma.
- Perforation of the aorta or renal artery wall
- Thromboembolism due to blood clotting or air bubbles

! This device should be exclusively applied by consultant cardiologists who are trained and experienced in vascular catheter interventions.

This training should be supervised and certified by an expert in the field of cardiovascular laser application.

Catheter Check

During cardiovascular applications, especially during catheterization of the renal or pulmonary artery and laser radiation, the **HypertenoLas®** is subjected to a variety of mechanical and thermal strains. If catheter damage is suspected, both the mechanical and optical integrity of the catheter must be verified by visual inspection.

Laser Safety

The **HypertenoLas®** must be used in compliance with current laser safety regulations. Operators and personnel must have sufficient knowledge about potential hazards and safety measures for medical laser application.

Safety Hazard

Sympathetic denervation is performed by means of invisible continuous wave laser light at 1064nm wavelength. Due to low absorption in water, it can easily penetrate the eyeball, and irreversible damage of the retina is the primary safety hazard.

!	Never point a connected <i>HypertenoLas</i>® towards a person or reflecting surfaces.
	Enable laser radiation only when the <i>HypertenoLas</i>® is positioned inside the body.
	Disabling of the laser footswitch is: 1. advisable after each laser application, is 2. mandatory before removing the catheter from the patient.
	Check the laser beam profile only by means of the pilot beam, never by releasing invisible 1064 nm laser light.
■	Use eye protection that blocks the 1064 nm wavelength.

Laser safety measures primarily include the prevention of inadvertent laser emission and eye protection for people staying inside the laser area (defined by the area, in which radiation exceeds the authorized value).

The *HypertenoLas*® causes a laser area of 1.3 m radius, measured from the catheter tip (calculated according to laser safety regulation DGUV-V 11 (previously BGV B2 for an irradiation at 20 W over 10 s).

Safety hazards are significantly reduced, if laser emission is inhibited on the laser panel whenever the *HypertenoLas*® is handled outside the patient. In this special case, a temporary laser area is caused inside the patient, when the foot switch is enabled. No hazardous radiation can be released into the operating room.

For laser application as described safety goggles are not necessary. However, laser goggles may be mandatory due to local safety regulations!

To visualize the course of the non-visible laser beam, medical lasers are equipped with a harmless red pilot beam. This beam is suitable to control the radiation characteristics of the *HypertenoLas*® and to study the different beam properties in air and water.

Medical lasers and accessories such as optical fibers or catheters are controlled by a series of hardware switches to avoid inadvertent radiation or radiation under inappropriate conditions.

Before laser application can be started by pressing the laser foot switch, the catheter must have been properly connected and the foot switch must have been enabled on the laser front panel. Make sure to regularly check the functions of the laser safety chain.

RISK MANAGEMENT

A *HypertenoLas*® that is inserted into the patient cannot radiate hazardous laser light into the environment. Laser hazard increases when the connected catheter is handled outside the patient, e.g., during visual inspection of the catheter. Therefore, it is recommendable to enable the laser footswitch only when the *HypertenoLas*® is already placed on the target site in the patient's body and to inhibit laser emission after laser application.

Hazards can be also provoked when the *HypertenoLas*® is handled inappropriately, especially when the optical cable is bent around sharp edges or is jammed under wheels or heavy weights.

The fiber tip of the *HypertenoLas*® mounted at a given distance from the end hole in the central lumen of the catheter is effectively protected against mechanical forces acting on the catheter during vascular catheterization.

However, it must be strictly avoided to clean its distal end hole by means of sharp objects such as metal pins or cannula, as this might cause fiber breakage. In this case, application of the *HypertenoLas*® would not be safe anymore, and it could endanger patients.

HANDLING

Preparation of the Laser System Calibration

The *HypertenoLas*® is optimized for radiation in an aqueous medium and cannot be calibrated in air. Light transmission is therefore measured by the manufacturer. Only qualitative control of the laser spot must be performed by operators (see chapter "Catheter Control").

Pump Triggering

Pump triggering via the laser footswitch should be checked before clinical use. Make sure that the footswitch changes the pump flow rate, and that laser light is only applied simultaneously with the increased flow from 15mL/min continuous flow to 35 mL/min lasing flow.

!	The <i>HypertenoLas</i>® must be irrigated and be surrounded by an aqueous medium. Laser radiation without irrigation can destroy the <i>HypertenoLas</i>® and could endanger patients.
---	--

PREPARATION OF THE PATIENT

The patients should be in a good clinical condition. Preferably, primary diagnosis is performed without anti-hypertensive medication. An additional venous access line is needed for volume and electrolyte substitution, antithrombotic and optional for other drug treatment.

Percutaneous venous or arterial puncture is performed under sterile conditions. In general, transluminal access to the renal arteries is achieved via a femoral artery. From there, the aorta abdominalis and the left and right renal arteries can be catheterized. For pulmonary denervation a femoral vein is punctured.

Connecting the Catheter

Main functions of the *HypertenoLas*® and its connections to external devices must be checked prior to a clinical application. This mainly includes:

- Integrity and sterility of the packaging,
- Catheter irrigation,
- Laser spot

Thoroughly inspect the packaging and expiry date before opening. Draw the cover sheet off the double sterilization envelope and remove the catheter. For connection of the optical cable, first remove the protective cap from the optical plug.

Be sure not to touch or soil the polished proximal end of the optical fiber. Insert the plug into the coupling port of the laser and fully tighten the screwcap until the safety interlock reacts by switching on the pilot laser.

Check the course of the optical cable. It should not run next to wheels or across sharp edges.

Avoid bending radii below 5 cm. Light should only be radiated from the catheter end hole. Do not use the catheter if light is leaking through the cable or the catheter body.

Plug the electrical cables in the ECG unit according to chapter "Connections" Insert sterile original tubing into the peristaltic pump.

Connect the inflow to heparinized saline and the outgoing line to the *HypertenoLas*®. Set the flow rate of the pump to

!	Inappropriate handling, bending around sharp edges, jammed under wheels or heavy weights can destroy the catheter.
---	---

35 mL/min and flush the line and catheter to remove air bubbles.


By releasing the footswitch, a continuous flow of 15 mL/min saline must leave the distal catheter end. Control the tightness of all connections. The tubing must not be distended by excessive pressure.

Now test the program! Set the continuous background flow rate to 15mL/min, and lasing flow rate to 35 ml/min. Start the background flow and step on the laser footswitch while laser emission is still inhibited.


The pump must now automatically switch to the lasing rate. When releasing the footswitch, the pump must return to the continuous background flow.

Eventually, inspect the laser spot. This test must be performed after irrigation, as the dry catheter radiates a completely different laser beam.

Dip the catheter tip into a white cup filled with sterile water or saline. Avoid reflective metal surfaces. The laser spot must be visible as a clear-cut continuous ring ("Donut"). An additional spot may be visible in its center that must not be brighter than the ring. Do not use catheters with anomalous laser spots.



The laser beam frontal view



The laser beam lateral view

Catheter Insertion

The following description of renal and pulmonary artery laser applications is a proposal based on general experience in the use of cardio-vascular catheters. Investigators may change procedure details according to their personal experience.

The *HypertenoLas*® is inserted percutaneously via a Guiding system, special long introducer sheaths and steerable guides such as the AGILIS (St. Jude). Guiding sheaths and introducers adapted to the *HypertenoLas*® are commercially available.

Do not use inappropriate introducer sheaths or guiding catheters as they may jeopardize the functions of the *HypertenoLas*® and increase the complication rate!

Procedure:

After aseptic skin preparation of the insertion area and sterile covering of the groin, a femoral vessel, a vein, or artery is punctured by using Seldinger technique.

A guide wire is introduced through the cannula into the - femoral artery and is advanced, under fluoroscopic guidance in the aorta abdominalis for perirenal or - via the femoral vein when pulmonary denervation is aimed at.

- If necessary, make a small incision to widen the puncture site for insertion of the tapered vessel dilator.

The proximal part of the guide wire must remain outside the patient over length of approx. 90 cm.

- The dilator in the guiding sheath is pushed over the wire and is advanced under fluoroscopic guidance into the abdominal aorta or pulmonary trunk respectively.
- Remove wire dilator while flushing with heparinized saline 0.9% 60mL/min via its sidearm.
- for RSD perform an aortogram via the guiding sheath, measure and calculate the vessel dimensions as described above and on page 4 and manipulate the guiding sheath from the aortic route towards the branch of the renal artery and position the distal end of the guide 5- 10mm from the renal artery ramification.

- Flush the **HypertenoLas®** fully to remove air bubbles and insert its tip into the hemostatic valve and push it into the guiding sheath.
- Advance the **HypertenoLas®** under continuous flushing with heparinized saline at 15mL/min through the guiding sheath and position the catheter tip under Xray control at the endhole of the guide.
- Start laser applications by using power of 15W, and time as calculated for 10J/mm²
- Deliver 2-3 application at intervals of 3-5s and repeat the procedure in the other renal artery, and if needed in - an accessory renal vessel too.
- Remove the **HypertenoLas®** and perform a final control aortogram.

For PSD after the pulmonary angiogram laser applications are performed as described, including LFS and HFS pace mapping. Pulmonary artery application times are limited to the effects achieved under HFS.

Catheter manipulation

The **HypertenoLas®** can be manipulated by advancing, withdrawing, and twisting the steerable guide and can be pushed up to 10 cm beyond the end hole of the Guide towards the targeted vessel segment.

Laser Application

The laser is used in the continuous wave irradiation mode at a distal laser power of 15 W. Energy setting must be applied according to the targeted inner surface of the vessel wall.

To avoid inadvertent release of laser light, triggering by the laser footswitch must be dis-enabled during catheter insertion and manipulation.

After Xray guided positioning of the catheter tip in the renal artery segment, laser radiation is enabled by means of the laser button on the laser front panel, and laser radiation can then be released by the two-staged laser foot switch.

By depressing the foot-switch the peristaltic pump flow rate increases to the preset lasing value. Simultaneously laser radiation is released. Laser alarm sounds during laser application.

Radiation is terminated automatically after the preset pulse duration or by releasing the footswitch. Longer application times are achieved by consecutive pulses. In this effort, the default pulse durations should be applied, and the radiation statistics displayed on the laser panel should be documented.

For safety reasons, the laser footswitch should be dis-enabled immediately after each laser application, no matter if the catheter remains in the patient for further treatments.

Disenable laser radiation before removing the catheter from the patient.


Laser radiation may produce a painful sensation of heat or pressure in the chest. Or pain in the kidney region. In rare cases analgesics may be required.

Catheter Control

Whenever suspecting a reduced catheter performance or a malfunction, the **HypertenoLas®** should be inspected. This includes its general appearance and the shape of the laser beam profile. Post-application laser beam inspection is analogous to the initial catheter check after unpacking.

Do not use catheters with deteriorated beam profile.

Radiation from a damaged optical fiber can cause severe complications such as vessel burning with crater formation, perforation, rupture of the vessel wall, and may produce thrombi with the risk of renal and systemic infarction with stroke.



Deteriorated beam profiles:

- 1. Interrupted ring (tip deformation casts a shadow): Reduced effect, potential overheating of the catheter**
- 2 Diffuse radiation (following contamination of the fibre tip): Reduced effect.**
- 3. Collimated and partly asymmetrical radiation (due to fibre breakage) DANGER! High risk of perforation; increased laser safety hazard.**

Connections

Laser

The **HypertenoLas®** is connected to a 1064 nm diode laser via a flexible optical fiber provided with a FSMA connector at its proximal end. The connector is provided with a chip as a time limiter of the laser transmission via the catheter. As light transmission and radiation characteristics of the catheter depend on the laser model, only the Laser:

CardioVasclLas® 1064, LPS
can relate to the **HypertenoLas®**

For protection of the fiber from inadvertent thermal damages caused by fiber contact with blood or tissue the Laser is provided with a Light-guide Protection System (LPS) that switches off the laser automatically in case of imminent overheating of the fiber tip.

There is no electrical connection between the laser and the patient or the working field. Damages of the Laser or the patient by a Cardioverter-Defibrillator can be ruled out.

Catheter Irrigation

The **HypertenoLas®** is flushed continuously with heparinized saline by means of a peristaltic pump. It is irrigated via a female Luer connector. Irrigation must be performed by a laser-triggered peristaltic pump. This is mandatory for a safe and correct application of the **HypertenoLas®**.

Irrigation is to be performed at a rate of at least 15 ml/min, rate that is automatically increased to 35ml/min during laser application. Manual pump adaptation is not recommendable. Do not use other pumps than the laser accessory peristaltic pump: **IriFlowLas®**

EXCLUSION OF LIABILITY

The **HypertenoLas®** is used in the thoroughly aggressive environment of the human body. Moreover, the needed flexibility and the very small diameter of the catheter invariably result in limited reliability.

Catheters can fail for any number of reasons, among which there are medical complications, rejection reactions of the body, fibrosis, dislocation, erosion or migrating through body tissue or resulting from a break or tear in insulation. Despite the greatest possible care taken in development, choice of components, assembly, and final control prior to delivery, the **HypertenoLas®** can become damaged by negligent handling or other influences, during or after introduction.

As a result, we do not assure or guarantee in any way that a disturbance or discontinuation of function will not occur, nor that the human body will not reject the introduced catheter, nor that there will be no medical complications, including myocardial perforation, because of introducing the catheter. Since the accessories can become damaged by negligent handling or other influences before, during or after introduction, we do not assure or guarantee in any way that a disturbance or discontinuation of function will not occur.

However, in cases of defects found prior to its use the no contaminated catheter can be replaced by the manufacturer, provided the catheter is sent back with documents describing the defect found.

The **HypertenoLas® is for single use only. Re-sterilization and reuse of the catheter would jeopardize its function and endanger the patient.**

The **HypertenoLas®** is sold in "as is" condition. No responsibility will be assumed for any deficiencies that are not immediately declared upon delivery of the goods.

The purchaser assumes the total risk related to the quality and function of the catheter and accessories when they are put into use.

LasCor® assumes no responsibility whatsoever for any loss, damage, or injury, be it directly or indirectly related to the catheter or accessories or determined to be subsequent damage resulting from the use thereof.

Consequently, **LasCor®** does not and will not assume any expenses incurred by the purchaser or a third party, ensuing from the use, malfunction or total failure of any catheter or accessories.

This exclusion encompasses physician's fees, costs of hospitalization, costs incurred by using medicinal products, any secondary expenses, and all subsequent damages.

Let it be known that no institution, organization, or individual has ever been empowered or in any other manner obtained the right to issue any notification deviating from the above or to make any guaranty in the name of **LasCor®**.

NOTE

After its use the catheter is biologically contaminated and must be disposed according to the local regulations. The catheter does not contain toxic components that would need special measures.

TECHNICAL DATA: **HypertenoLas®**

Body size	8 F
Useable length	115 cm
Total Length	3 m
Optical fiber Core	400 µm
Numerical aperture	0.22
Beam divergence fiber in water	70°
Distributor	3way
Irrigation tube	15 cm
Laser-connector	FSMA
Catalogue No.:	H 003-115-XXX
UDI-DI:	4260691560030
Packaging:	carton box with one sterile set
Delivery box:	with 10 single carton boxes


SYMBOLS

This medical product is labelled according to **DIN EN ISO 15223-1:2022-02** by using the following graphical symbols:


	Medical Device
	Manufacturer
	Lot number
	Date of manufacturing
	Sterile, sterilization method
	For single use only
	Regard operating manual
	To be used until

Peri-renal artery nerve modulation

Example for calculation of energy density to be applied for peri-renal nerve modulation in the treatment of arterial systemic hypertension.



Left renal angiogram **showing** an **18.2mm** long vessel length from its ramification, with a uniform caliber of the renal artery and a diameter of **D = 4.2mm**.



Calculated **CIRCUMFERENCE**, diameter multiplied with 3.14 (π): $5 \times 3.14 = 15.7\text{mm}$.

Based on experimental results: *J Vet Sci Ani Husb* 9(1):103-113 2021 ISSN: 2348-9790 for optimal peri-renal nerve modulation, for a safe and effective permanent perivascular denervation (RSD) laser applications at $10\text{J}/\text{mm}^2$ is aimed at.

Calculation:

- The **diameter** is 5mm, the **circumference** is $5 \times 3.14 (\pi) = 15.7\text{mm}$.
- The **length** of the tubular illuminated vessel inner wall is approximately 5.0mm
- The illuminated sheath surface is circumference x length (5mm): $15.7 \times 5 = 78.5 \text{mm}^2$.
- Laser application at **15W** on 78.5mm^2 will achieve: $15:78.5 = 0.19 (0.20) \text{W}/\text{mm}^2$.

To achieve $10\text{J}/\text{mm}^2$ an irradiation time of 50s is needed: $0.19 \text{W}/\text{mm}^2 \times 50\text{s} = 9.51\text{J}/\text{mm}^2$
(Approximated $0.20 \text{W}/\text{mm}^2 \times 50\text{s} = 10\text{J}/\text{mm}^2$)

To achieve approximately $10\text{J}/\text{mm}^2$ you must adapt radiation times to the vessel diameter you have measured on your screen. See table below.

Artery Diameter (mm)	Radiation Time (s)	
	Calculated	Suggested
2,6	27	2 x 15
3,0	31	2 x 15
3,5	37	2 x 15 + 10
4,0	41	2 x 15 + 10
4,5	47	3 x 15
5,0	52	3 x 15 + 5
5,5	58	4 x 15
6,0	63	4 x 15
6,5	66	4 x 15

Procedure


Laser power is preset to **15W**. Adopted/assumed length of illuminated sheath is **5mm**. To achieve $10\text{J}/\text{mm}^2$ in a 3.0mm diameter vessel laser application should start with **2x15s** with an interval of 3s.

If longer radiation times are needed to achieve the aim of $10\text{J}/\text{mm}^2$, e.g., vessel diameter is 3.5 or 4.0mm, three seconds after the second 15s pulse **10s of radiation** must be added for an effective renal perivascular nerve modulation.

There is a stepwise increase of radiation time and a linear increase of the level of total energy applied.


Peri-pulmonary artery nerve modulation

Laser effects can be controlled by laser catheter LFS and HFS stimulation. The level of energy to be applied for pulmonary artery perivascular nerve modulation in the treatment of resistant pulmonary artery hypertension depends on the results of LFS and HFS.



ECG recordings during PA stimulation mapping.

- No atrial nor ventricular capture during LFS at a point in the right PA. HFS results in significant slowing of atrial rhythm. Cycle length prolongation from 617 to 1027ms.
- The ablation catheter is positioned at the same point as in the A tracing. After energy application no change in the atrial cycle length during HFS.



Note: Catheter orientation is not a major determinant for lesion formation. In flat catheter position even substantially larger sizes of lesion can be achieved by application of longer radiation times.

ANNEXES:

- 1 Patient Information
- 2 Patient's statement written informed consent

The patient information sheet must be handed out to and must be signed by the patient together with the informed written consent.

- 3 Fax Template: Post Market Clinical Follow-up is compelling after the use of the catheter.

Version 2023-03-26

LasCor GmbH, Laser – Medical Devices
No 4, Schlesierstraße
D-82024 Taufkirchen, Germany

Tel.: +49 (0)89 759 5596
Fax: +49 (0)89 759 5770

E-Mail: info@lascor.de www.lascor.de (With IFUe)

Patient Information

Renal and Pulmonary Laser Denervation (LD)

What it's all about?

Introduction

Perivascular Laser Denervation (LD) is a non-surgical technique, a percutaneous, transluminal, minimally invasive catheter procedure that selectively inactivates perivascular sympathetic nerves causing your high blood pressure.

It involves percutaneous puncture or a tiny incision in the skin and is done without putting you to sleep. In fact, you will be awake during the procedure and be able to talk with the staff. The procedure is used for treating resistant systemic or pulmonary hypertension that has not responded adequately to medication and your improved lifestyle. Catheterization of the renal or pulmonary arteries with the catheter enables doctors to destroy or modify selectively the perivascular nerves by laser radiation.

Current catheter ablation techniques utilize radiofrequency (RF) energy or ultrasound (US). By using the open-irrigated catheter **HypertenoLas** laser ablation is a new promising method developed with the intention to reduce the risks, to increase success rates, and to shorten procedure duration times for renal or pulmonary sympathetic denervation.

For the laser procedure, doctors insert after a small incision of a vessel, or percutaneously, a steerable guiding sheath for vessel angiography to visualize vessel anatomy and localize its root for selective catheterization. After angiography the laser catheter is introduced through that steerable sheath and is advanced under X-ray control together with the sheath into the targeted vessel, the main renal artery up to its ramifications or the pulmonary trunk, and then pass laser light through the catheter. The irradiated circular area of the vessel at lengths of several millimetres is repeated if needed. If the treatment is successful, you will be cured from your hypertension or ineffective medication will work. The laser application is painless; however, the method is experimental but has shown to selectively inactivate the sympathetic nerves without damaging renal or pulmonary artery wall.

In case if you suffer also from cardiac arrhythmias that were ablated prior to renal denervation, renal as well as pulmonary sympathetic laser denervation can be performed also by using the catheter **RytmoLas**® in the same session after the arrhythmia ablation (Patient Information "Laser Catheter ablation of Cardiac arrhythmias").

Alternate Therapies to Renal or Pulmonary Sympathetic Laser denervation (SLD)

The following antihypertensive therapies are nowadays available:

1. *Medication* (which was not successful in the treatment of your resistant hypertension).
2. *Baroreceptor stimulation* (is a surgical procedure)
3. *Radiofrequency and Ultrasound Catheter Ablation* (also experimental procedures).

Preparing for the Ablation Procedure

Unless you are already hospitalized, you will probably be admitted to the hospital or, in some cases, you may undergo the treatment in outpatient hypertension unit.

Several routine laboratory tests will be performed prior to the intervention including an ECG, blood tests, and possibly ambulatory long-term blood pressure measurements. These may be done one or two days ahead of catheter ablation. The doctor performing the ablation procedure will review your medical history and examine you. You may be seen by the doctor at the office several days before the procedure.

The doctor will explain the ablation technique, its purpose, potential benefits, and possible risks. This is a good time to ask questions and, most importantly, to share any feelings or concerns you may have about the ablation intervention.

You will be asked to sign a consent form, a document that allows drugs to be injected directly into the vein if necessary. To help you relax, you will be given a sedative if necessary. Generally, you will be asked not to eat 6 hours or drink two hours anything prior to the intervention but if you are scheduled for an afternoon study, you may be given a light breakfast. You may have a small amount of water with medication.

You may be asked to stop taking certain medications for two or three days before the ablation procedure. Be sure to check with your doctor several days before the intervention. Bring a list of all the medications you are currently taking. It is very important for the doctor to know the exact names and dosages of any medication that you take. Be sure to mention to the doctor if you have had allergic reactions to any medication.

For your comfort, empty your bladder as completely as possible before the procedure starts. Once preparations are completed, you will be taken to the catheter laboratory where the procedure will be performed. You will be transferred to an X-ray or NMR table. The table has a large camera above it and television screens close by. The equipment in the lab also includes heart monitors and various instruments and devices. The lab team generally includes the angiologist with special training, radiologist, an anaesthetist, an assistant, nurses, and technicians.

During the Ablation Procedure

After being positioned on the table, you will be hooked to a variety of monitors and manifolds, and you will be covered with sterile sheets. The staff will be wearing sterile gowns and gloves. The groin where the catheter will be inserted is cleansed thoroughly. A local anaesthetic is injected into the skin with a tiny needle to numb the area. This may cause a stinging sensation. A small incision is made in the skin, and a needle system is used to puncture the artery, into which the catheter will be inserted. The open-irrigated laser ablation catheter **HypertenoLas®** is a long flexible tube provided with a flexible optical fibre mounted coaxially in its inner lumen that can transmit laser light. The catheter is not heated up and it is flushed continuously with saline through its central lumen during the procedure.

The catheter will be advanced into the main renal artery or pulmonary artery depending on your type of hypertension, systemic or pulmonary to irradiate this artery at a length of approximately 5mm. Laser applications at 15W can be repeated depending on the inner surface of the targeted vessel wall. By doing so perivascular sympathetic nerves are selectively inactivated whereas the vessel wall and perivascular tissues are less absorbing laser photons and therefore are not damaged permanently by the laser induced heat. This process of reducing sympathetic activity of periarterial nerves will reduce your blood pressure, in optimal case to normal values permanently and without the need of any antihypertensive drugs.

You will be awake during the procedure; although medication will be given to help you relax. The staff will be monitoring your process constantly. Let the staff know at any time if you experience pain or discomfort. The laser procedure usually is not painful, although you may feel some discomfort during the insertion of the catheters or during laser application.

There may also be some discomfort from lying still for a long time. You will not feel the catheters moving through the blood vessels, the aorta, or pulmonary arteries. Depending on the anatomy of your arteries procedures may last from one to two hours.

The insertion of catheters is accompanied by certain risks. Some patients may develop bleeding at the insertion site. Blood collects under the skin resulting in local swelling and "bruise". Both swelling and bruise will disappear in time as the blood is slowly absorbed by the body. Less frequently, ablation procedures may be associated with more serious complications.

These include damage to blood vessels, formation of blood clots and infection. Fatalities such as renal lesions of coagulation necrosis cannot occur provided laser applications are not performed in small artery ramifications close to the kidney hilum. However, laser procedure-related complications cannot be completely ruled out and of the risks. To learn about your particular risk, you should discuss the matter with the doctor.

After the Ablation Procedure

After the laser procedure is completed and the catheters are removed, the doctor or nurse will apply slight pressure over the groin for about 10 to 20 minutes. This is done to prevent bleeding. The doctor may close the incision with a few stitches.

Subsequently, a dressing with pressure will be applied over the wound or puncture. You will then be transported to your room or to the recovery area. You will probably be allowed to drink and eat following the procedure but check with the nurse. You will lie flat in bed for six to eight hours to allow a small seal to form over the puncture in the artery.

During this time, you may move your foot or wiggle your toes. You may move your arms freely. You will stay under observation for one night and if you are able to go home after control of the punctured groin or wound. However, if the procedure was without complications your doctor may decide to perform the treatment ambulatory.

Please remember the following to ensure a quick recovery:

- Limit your activities during the first 24 hours at home. You can move about, but do not strain or lift heavy objects.
- If you notice new blood on the dressing, place your fingers over the site and press for about 20 minutes. If bleeding continues, call your doctor, or the ambulance to the nearest emergency room while continuing to apply pressure.
- You may remove your dressing three days after the procedure and have shower but wait another day or two before taking a bath. You
- A black-and-blue mark (bruise) or a small lump under the skin at the insertion site is common. These generally disappear within three to four weeks.
- Call your doctor if the insertion site becomes painful or warm, or if it shows any other signs of infection.
- Ask your doctor when you can return to normal activities, usually one week after procedure, and if there are any specific restrictions.
- Ask your doctor about your medication, which ones to continue and which ones to stop.

Follow-up Ablation Procedure

If successful, the laser treatment has inactivated a part of the perivascular renal or pulmonary nerves, thereby reducing your blood pressure. Hopefully you will be cured from hypertension permanently without the need for medication. In rare cases, medication may still be needed after the ablation procedure, and medication not effective prior to the laser treatment may now work. It can be anticipated that your blood pressure will further decrease in the months following the laser treatment.


Follow-up examinations can be performed in the outpatient clinic. If results indicate that medication is still required, needs adjustment or is not effective, you may be brought back to the laboratory to repeat the arrhythmia ablation attempt.

For your *long-term follow-up*, you will be seen by your doctor in his office or at the outpatient department of the hospital for a regular control including physical examination and blood pressure control, preferably long-term, at time intervals decided by your doctor.


However, you should come to your doctor whenever you have other signs or symptoms which in your opinion may be related to the hypertension, as you have suffered from prior to the ablation procedure. Eventually, regular controls of your blood pressure will be performed by a visit at your house doctor.

Don't hesitate and ask your doctor if you need more information or if some of the above is still unclear to you or if you would like to know more about the procedure, possible complications, and risks.

Your doctor will give you further information.



Anterior posterior Xray image showing the position of the laser catheter in the pulmonary artery trunk during laser application (arrow). The lower arrow points towards the guiding Catheter in the inferior vena cava.



CATHETER ACCESS IN THE RIGHT GROIN

Scheme showing the main renal artery segment where laser catheter radiation is aimed at, and the artery is denervated.

Patient's Statement and Written Informed Consent

Signature File in accord with 21 CFR 50 with the Patients Statement

I had detailed discussion(s) with my doctor (print): _____
Concerning my questions, problems, concerns and doubts of the laser ablation procedure and the possible complications and risks, I have received a complete file of **Patient Information** (four pages) and I have no further questions.

Herewith I agree with the proposed Renal Pulmonary Laser catheter denervation, with the emergency interventions, necessary for the treatment of possible complications, and with the follow-up controls or studies.

I do not give my consent for the laser catheter ablation because (comment):

.....
.....
.....
.....
.....



Delete what not applicable.

Patient: _____
Name (print)

Date

Signature

Investigator: _____
Name (print)

Date

Signature

